SURG-I-BAND®

Color Coding

Data Matrix Barcode

MedEnvoy Global B.V.
Prinses Margrietplantsoen 33
Suite 123
2595 AM The Hague
The Netherlands

SCANLAN ★

INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921

SCANLAN INTERNATIONAL, INC.
ONE SCANLAN PLAZA
SAINT PAUL, MINNESOTA 55107 U.S.A.
800-328-9458 • 651-298-0997
FAX 651-298-0018
www.scanlaninternational.com

Surg-I-Band® Color Coding Data Matrix Barcode

Clinical Benefits:

Surg-I-Band® color coding is applied to surgical instruments and devices for identification purposes. Surg-I-Band® is made of a flexible polymer that can be used on all shapes and sizes of surgical instruments. This product may be used with steam, flash, and ethylene oxide (EO) sterilization.

Surg-I-Band® Data Matrix is offered in various colors and sizes. They are sequentially numbered with 2D barcodes and conveniently packaged on a spool. For added ease of use, the Surg-I-Band® Dispenser, Carousel, and Remover are optional accessories sold separately.

Intended User:

The intended users of Surg-I-Band® are Sterile Processing Department Personnel and associated areas within the health care facility setting.

Intended Use:

Surg-I-Band® color coding is intended for coding and identification of surgical instruments.

Instructions for Use:

1. Surgical instruments and hands must be dry and free of residue.
2. Select the desired color.
3. Keeping the data matrix in the center, disperse enough Surg-I-Band® to overlap itself once but not more than twice around the selected site on the surgical instrument.
4.  **CAUTION:** Do not overlap the data matrix barcode.
5. Place data matrix barcode in a visible location on a flat surface (minimum shaft diameter of 1/8").
6. Surg-I-Band® data matrix on spools must be cut with a sharp edge. The Surg-I-Band® Dispenser is recommended for the best application.
7.  **CAUTION:** Surg-I-Band® must not be stretched (stretching will weaken Surg-I-Band® bonding capability).
8.  **CAUTION:** Place Surg-I-Band® data matrix on the finger rings of scissors. It is not recommended that Surg-I-Band® data matrix be applied to the shanks of scissors, as this may interfere with the operation of the surgical instrument.
9. Color combinations may be used on the same surgical instrument. Place additional color coding so it does not block the data matrix barcode or the white border surrounding the barcode.
10. Upon completing the application of Surg-I-Band® color coding, autoclave the surgical instrument(s). Surg-I-Band® color coding will bond itself to the surgical instrument(s).
11. Inspect Surg-I-Band® and remove if it shows wear.
12. Use the Surg-I-Band® Remover to assist in the removal of Surg-I-Band® and clean with facility approved surgical instrument cleaner.
13. Use the Surg-I-Band® Carousel for organization and convenient storage.

Contraindications:

Surg-I-Band® is not intended for implantation.

Surg-I-Band® is not intended for re-use if damage is observed.

Surg-I-Band® is not intended for use in an EO cycle until after adhered to substrate by autoclave cycle.


Surg-I-Band® is not intended for use in dry heat sterilization.

Surg-I-Band® is not intended to be placed where it will interfere with device function.

Warnings/Cautions:

 **CAUTION:** Surg-I-Band® is not intended for use in dry heat sterilization.

 **CAUTION:** Surg-I-Band® is not intended for use in hydrogen peroxide gas plasma.

 **CAUTION:** Destroy Surg-I-Band® that might have been contaminated with prions (Creutzfeldt-Jakob Disease, CJD). Do not reuse.

Risk of Reuse Statement:

Reprocessing or re-use of Surg-I-Band® is permitted if this IFU is followed; requiring thorough inspection each time a surgical instrument or device is reprocessed and removed, if necessary. Any misuse of Surg-I-Band® shall void any warranty or liability by the manufacturer for injury or loss associated with inappropriate use.

Requirements for Special Facilities, Training, or Particular Qualifications of the User and/or Other Persons:

Not Applicable

Methods for Safe Disposal:

The user is responsible to dispose of the device and packaging according to local and national requirements.

Incident Reporting Notice:

Any serious incident that has occurred in relation to the device should be reported to Scanlan International, Inc. and the competent authority of the Member State in which the user and/or patient is established.

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