


A/C LOCATOR® U-Shaped Graft Markers
RADIOMARK® Disk-Shaped Graft Markers



INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921

Designed For Assessment of Graft Patency During Postoperative Angiography



A/C Locator® U-Shaped Graft Markers

The A/C Locator® U-Shaped graft marker is designed to be applied with 5-0 or 6-0 suture after placement of coronary and peripheral grafts.

Manufactured with lightweight medical grade silicone, it will not produce a reflective glare when used with a CT scan or MRI. It offers flexibility, biocompatibility, meets all implant criteria and will not migrate. Radiopaque.

Catalog Number
1001-83



(shown actual size)

Radiomark® Disk-Shaped Graft Markers


The Radiomark® disk-shaped graft marker is also designed to be applied with 5-0 or 6-0 suture after placement of coronary and peripheral grafts.

Manufactured with lightweight medical grade silicone, it will not produce a reflective glare when used with a CT scan or MRI. It offers flexibility, biocompatibility, meets all implant criteria and will not migrate. Radiopaque.

Catalog Number
1001-831



(shown actual size)



Features Include...

- Unique design
- Lightweight
- Biocompatible
- Medical grade silicone
- Radiopaque
- Not made with natural latex rubber
- Tissue compatible
- Meets all implant criteria

For more information or to place an order,
please contact your local Scanlan representative,
or call... 800-328-9458
FAX: 651-298-0018
www.scanlaninternational.com



SCANLAN ★

INTERNATIONAL