

Extracranial-Intracranial (EC/IC) Bypass Instrumentation



SCANLAN ★

INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921
800-328-9458 • 651-298-0997 • FAX 651-298-0018 • www.scanlangroup.com

USA



INTERNATIONAL

Surgical products designed and manufactured by the Scanlan family since 1921

EC / IC Bypass Instrumentation Set

Three lengths of instruments have been developed to accommodate the requirements for the varying depths of cranial bypass.


Additional features of this offering include:

- A precise 0.3 mm tip
- Scanlan's exclusive Diamond Dust™ technology to control needle rotation and enable atraumatic handling of needles, suture and tissue
- Streamlined tips on all forceps and needle holders provide enhanced visibility of micro needle and sutures
- Titanium forceps eliminate needle magnetism
- Metal instrument trays prevent damage to instruments during transport, cleaning and sterilization

SCANLAN®
EC / IC Bypass Instrument Set 9999-166

Total set includes three subsets
9999-162 (Short) pages 4 - 5
9999-163 (Medium) pages 6 - 7
9999-164 (Long) pages 8 - 9

EC / IC Bypass Instrument Set may be purchased as a set or as individual components.



0.3 mm 45° Tying Forceps



0.3 mm Gerald Type Taper Straight Forceps



0.3 mm Gerald Type Taper Curved Forceps



1 mm Ring Tip Forceps



0.3 mm Curved Needle Holder



45° Vannas Micro Scissors



SCANLAN®
EC / IC Bypass Instrumentation Set
 Short Instrument Set 9999-162

SCANLAN®
Forceps
 Round Handle



Instrument is shown approximate size

SCANLAN® Micro Suture Tying Forceps, Titanium

Qty.	Cat. No.	Description	Total Length
1	3003-154-3	0.3 mm Straight Platform Gerald Type Taper DIAMOND DUST™	6" (15 cm)



Instrument is shown approximate size

SCANLAN® Micro Suture Tying Forceps, Titanium

Qty.	Cat. No.	Description	Total Length
1	3003-156-3	0.3 mm Curved Platform Gerald Type Taper DIAMOND DUST™	6" (15 cm)



Instrument is shown approximate size

SCANLAN® Micro Forceps, Titanium

Qty.	Cat. No.	Description	Total Length
1	3003-164	1 mm Straight DIAMOND DUST™ Ring Tip	6" (15 cm)



Instrument is shown approximate size

SCANLAN® Micro Suture Tying Forceps, Stainless Steel

Qty.	Cat. No.	Description	Total Length
1	4004-03-3A	0.3 mm Angled 45° Angle DIAMOND DUST™ Tips	5 1/2" (14 cm)

SCANLAN®
EC / IC Bypass Instrumentation Set
 Short Instrument Set 9999-162

SCANLAN®
Needle Holders
 Round Handle



Instrument is shown approximate size

SCANLAN® Castroviejo Micro Needle Holder, Stainless Steel

Qty.	Cat. No.	Description	Recommended Suture Size	Total Length
1	6006-56-3	0.3 mm Tapered Curved DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet	8-0 and Smaller	5 3/4" (14.5 cm)

SCANLAN®
Jacobson Vannas Micro Scissors
 Flat Handle



Instrument is shown approximate size

SCANLAN® Jacobson Vannas Micro Scissors, Stainless Steel

Qty.	Cat. No.	Description	Total Length
1	7007-140	45° Angle Micro Fine Blades	7" (18 cm)



SCANLAN®
EC / IC Bypass Instrumentation Set
 Medium Instrument Set 9999-163

SCANLAN®
Forceps
 Round Handle



Instrument is shown approximate size

SCANLAN® Micro Forceps, Titanium

Qty.	Cat. No.	Description	Total Length
1	3003-158-3	0.3 mm Straight DIAMOND DUST™ Gerald Type Taper Platform Tips	7 1/4" (18.5 cm)



Instrument is shown approximate size

SCANLAN® Micro Forceps, Titanium

Qty.	Cat. No.	Description	Total Length
1	3003-159-3	0.3 mm Curved DIAMOND DUST™ Gerald Type Taper Platform Tips	7 1/4" (18.5 cm)



Instrument is shown approximate size

SCANLAN® Micro Forceps, Titanium

Qty.	Cat. No.	Description	Total Length
1	3003-164	1 mm Straight DIAMOND DUST™ Ring Tips	6" (15 cm)



Instrument is shown approximate size

SCANLAN® Micro Forceps, Stainless Steel

Qty.	Cat. No.	Description	Total Length
1	4004-97-3A	0.3 mm Angled DIAMOND DUST™ 45° Angle	5 1/2" (14 cm)

SCANLAN®
EC / IC Bypass Instrumentation Set
 Medium Instrument Set 9999-163

SCANLAN®
Needle Holders
 Round Handle



Instrument is shown approximate size

SCANLAN® Jacobson Micro Needle Holder, Titanium

Qty.	Cat. No.	Description	Recommended Suture Size	Total Length
1	3003-177-3	0.3 mm Tapered Curved DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet	8-0 and Smaller	7 1/4" (18.5 cm)

SCANLAN®
Jacobson Vannas Micro Scissors
 Flat Handle



Instrument is shown approximate size

SCANLAN® Jacobson Vannas Micro Scissors, Stainless Steel

Qty.	Cat. No.	Description	Total Length
1	7007-140	45° Angle Micro Fine Blades	7" (18 cm)

See pages 10-11 for optional SCANLAN® V.I.P.™ Sterilization Trays



SCANLAN®
EC / IC Bypass Instrumentation Set
 Long Instrument Set 9999-164

SCANLAN®
Forceps
 Bayonet Style | Flat Handle



SCANLAN® Micro Forceps, Titanium

Qty.	Cat. No.	Description	Total Length
1	3003-842	0.3 mm Straight DIAMOND DUST™	8 3/4" (22 cm)



SCANLAN® Micro Forceps, Titanium

Qty.	Cat. No.	Description	Total Length
1	3003-841	0.3 mm 45° Angle DIAMOND DUST™	8 3/4" (22 cm)



SCANLAN® Micro Forceps, Titanium

Qty.	Cat. No.	Description	Total Length
1	3003-840	0.3 mm Curved Up DIAMOND DUST™	8 3/4" (22 cm)



SCANLAN® Micro Forceps, Stainless Steel

Qty.	Cat. No.	Description	Total Length
1	3003-851	1 mm Straight DIAMOND DUST™ Ring Tip	8 3/4" (22 cm)



SCANLAN®
EC / IC Bypass Instrumentation Set
 Long Instrument Set 9999-164

SCANLAN®
Micro Needle Holders
 Bayonet Style | Round Handle



SCANLAN® Micro Needle Holder, Stainless Steel

Qty.	Cat. No.	Description	Recommended Suture Size	Total Length
1	3789-14-3	0.3 mm Curved Up DIAMOND DUST™ Jaws Regular Box Lock Without Ratchet	8-0 and Smaller	9" (23 cm)

SCANLAN®
Micro Scissors
 Bayonet Style | Flat Handle



SCANLAN® Micro Scissors, Stainless Steel

Qty.	Cat. No.	Description	Total Length
1	3788-02	Vannas 45° Angled Up On Flat Micro Fine Blades	8 3/4" (22 cm)

See pages 10-11 for optional SCANLAN® V.I.P.™ Sterilization Trays



Optional
SCANLAN® V.I.P.™ Sterilization Trays



2138-03

Description

Silver-gray color | Scan-Guard™ Silicone Inserts
 10" x 10" x 1.5" (25.4 x 25.4 x 3.8 cm)
 Capacity: 8-10

2140-03

Description

Silver-gray color | Scan-Guard™ Silicone Inserts
 10" x 15" x 1.5" (25.4 x 38 x 3.8 cm)
 Capacity: 10-12



Optional
SCANLAN® V.I.P.™ Sterilization Trays



2141-03

Description

Silver-gray color | base with accessory
 compartment | Scan-Guard™ Silicone Inserts
 10" x 15" x 2.6" (25.4 x 38 x 6.6 cm)
 Capacity: 15+



EC / IC Bypass Instrumentation Set



For more information or to place an order,
please contact your local Scanlan representative, or call...

U.S. & Canada: 800-328-9458

International: 651-298-0997

FAX: 651-298-0018

www.scanlaninternational.com



Surgical products designed and manufactured by the Scanlan family since 1921
800-328-9458 • 651-298-0997 • FAX 651-298-0018 • www.scanlangroup.com

USA



MK-B-ECIC Rev. B