

SCANLAN®

# Scanturian® Minimally Invasive Surgical Instruments

Precision | Performance | Durability | Quality

Individually handcrafted with meticulous detail to assure superior quality, precision and performance!

## — NEEDLE HOLDERS —

Carbon-Inlay Jaws



## — NEEDLE HOLDERS —

DIAMOND DUST™ Jaws



## — FORCEPS —

DeBakey Jaws | Resano Jaws



## — SCISSORS —

Straight | Curved | Angled



Highest Quality Surgical Instruments  
Since 1921



Joseph Scanlan Adam Scanlan Tim Scanlan Russell Scanlan Chenoa Scanlan Ron Peick

Dear Customer,

Since 1921, Scanlan International has been dedicated to working closely with surgeons to provide the highest quality, innovative surgical instruments.

Over the years, the specialties of cardiac, thoracic and vascular surgery are continuing to evolve from a conventional open approach to a minimally invasive surgical approach. This evolution to minimally invasive surgery (MIS) has been the catalyst for the next generation of specialty surgical instrumentation.

Scanlan continues to be a leader in the design and development of MIS specialty instrumentation. The new Scanturian® minimally invasive surgical instruments are individually hand-crafted to assure a superior degree of instrument quality, control and detail you expect in your Scanlan instruments.


Thank you for your continued support and demand for the finest surgical instrumentation. If we can ever be of assistance to you, please contact us.

Sincerely,

Timothy M. Scanlan  
President and Chief Executive Officer

Russell T. Scanlan  
Chief Operating Officer

# MIS SCANLAN® Scanturian® Minimally Invasive Surgical Instruments



### ScanShield® Shaft Feature

The ScanShield® shaft feature surrounds the hinge mechanisms to avoid suture entrapment and intraoperative trauma to valve tissue and chordae tendineae

### Tapered Shaft

The tapered shaft design allows for more robust moving mechanisms and efficient cleaning due to the increased shaft diameter at the distal tip

### Scanturian® MIS Scissors with SUPER CUT™ Blades

A unique blade design - combining a razor edge with a scissors edge, these scissors produce a clean, forward cut, eliminating pinching or slipping of tissue

### Instrument Identification

Color-coded and etching for easy intraoperative instrument identification



### Titanium Handle

Well-balanced, durable and handcrafted MIS Instrument design

### Handle Knurling

Proven handle knurling design provides comfortable, secure grip and ergonomic handle design


### Proximal End Cleaning Port

Luer lock flush port system allows for reliable automated and intraoperative instrument flushing

ScanShield® Patent no. 11,478,265



MIS **SCANLAN®**  
**Scanturian® MIS Forceps**



Scanturian®  
**MIS Forceps**

- DeBakey and Resano Jaws
- Well-balanced instrument
- High-quality and durable
- Securely grasps tissue
- Unique feature - SCANLAN® ScanShield® shaft

MIS **SCANLAN®**  
**Scanturian® MIS Forceps**

Cat. No.	Description	Distal Tip	Working Length	Total Length	Handle Color
9019-420-17	SCANLAN® Scanturian® MIS DeBakey Forceps 1.5 mm Straight Jaws Without Ratchets		6-3/4" (17 cm)	12-1/4" (31 cm)	
9019-420-25			9-7/8" (25 cm)	15-1/3" (39 cm)	
9019-430-17	SCANLAN® Scanturian® MIS DeBakey Forceps 2.5 mm Straight Jaws Without Ratchets		6-3/4" (17 cm)	12-1/4" (31 cm)	
9019-430-25			9-7/8" (25 cm)	15-1/3" (39 cm)	
9019-440-17	SCANLAN® Scanturian® MIS Resano Forceps 2.5 mm Straight Jaws Without Ratchets		6-3/4" (17 cm)	12-1/4" (31 cm)	
9019-440-25			9-7/8" (25 cm)	15-1/3" (39 cm)	

ScanShield® Patent no. 11,478,265

Please Note...  
Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request. For a complete listing of models visit... [www.scanlaninternational.com](http://www.scanlaninternational.com)






# SCANLAN® Scanturian® MIS Needle Holders



9019-610-17  
Actual Size



9019-610-25  
Actual Size



## Scanturian® MIS Needle Holders

- Curved and straight carbon-inlay jaws
- Well-balanced instrument
- High-quality and durable
- Needle holder securely holds the needle
- Tapered shaft
- Recommended for 2-0 and 4-0 suture size and smaller

## MIS SCANLAN® Scanturian® MIS Needle Holders

Cat. No.	Description	Distal Tip	Suture Size	Working Length	Total Length	Handle Color
9019-610-17	SCANLAN® Scanturian® MIS Needle Holder Curved Carbon-Inlay Jaws with Ratchet		2-0 and Smaller	6-3/4" (17 cm)	12-1/4" (31 cm)	
9019-610-25				9-7/8" (25 cm)	15-1/3" (39 cm)	
9019-614-17	SCANLAN® Scanturian® MIS Needle Holder Straight Carbon-Inlay Jaws with Ratchet		4-0 and Smaller	6-3/4" (17 cm)	12-1/4" (31 cm)	
9019-614-25				9-7/8" (25 cm)	15-1/3" (39 cm)	
9019-630-17	SCANLAN® Scanturian® MIS Needle Holder Straight Carbon-Inlay Jaws with Ratchet		4-0 and Smaller	6-3/4" (17 cm)	12-1/4" (31 cm)	
9019-630-25				9-7/8" (25 cm)	15-1/3" (39 cm)	

Please Note...  
Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request. For a complete listing of models visit... [www.scanlaninternational.com](http://www.scanlaninternational.com)






# SCANLAN® Scanturian® MIS Needle Holders



9019-640-17  
Actual Size



9019-640-25  
Actual Size



## Scanturian® MIS Needle Holders

- Straight DIAMOND DUST™ jaws
- Well-balanced instrument
- High-quality and durable
- Needle Holder securely holds the needle
- Recommended for 5-0 suture size and smaller

## SCANLAN® Scanturian® MIS Needle Holders



Cat. No.	Description	Distal Tip	Suture Size	Working Length	Total Length	Handle Color
9019-640-17	SCANLAN® Scanturian® MIS Needle Holder Straight DIAMOND DUST™ Jaws with Ratchet		5-0 and Smaller	6-3/4" (17 cm)	12-1/4" (31 cm)	
9019-640-25				9-7/8" (25 cm)	15-1/3" (39 cm)	

MIS Needle Holders with Scanlan® Diamond Dust™, a proprietary coating for needle holders provides a secure grip to prevent needle rotation and burring. Diamond Dust™ surfaces are guaranteed for 12 months and will be replaced at no charge if it wears away prior to that time.

Please Note...  
Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request. For a complete listing of models visit... [www.scanlaninternational.com](http://www.scanlaninternational.com)






# SCANLAN® Scanturian® MIS Scissors



9019-710-17SC  
Actual Size

9019-710-25SC  
Actual Size



## Scanturian® MIS Scissors

- Well-balanced instrument
- High-quality and durable
- Super Cut™ Scissors cut cleanly and without slipping while the handle provides tactile feedback
- SCANLAN® ScanShield® shaft

ScanShield® Patent no. 11,478,265

**Please Note...**  
Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible. Reflecting this unique process, all stated dimensions are approximate.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.  
For a complete listing of models visit... [www.scanlaninternational.com](http://www.scanlaninternational.com)



## SCANLAN® Scanturian® MIS Scissors

Cat. No.	Description	Distal Tip	Working Length	Total Length	Handle Color
9019-710-17SC	SCANLAN® Scanturian® MIS Super Cut™ Scissors Straight Fine Blades, 1 cm, Blunt Tips		6-3/4" (17 cm)	12-1/4" (31 cm)	
9019-710-25SC			9-7/8" (25 cm)	15-1/3" (39 cm)	
9019-715-17SC	SCANLAN® Scanturian® MIS Super Cut™ Scissors Curved Fine Blades, 1 cm, Blunt Tips		6-3/4" (17 cm)	12-1/4" (31 cm)	
9019-715-25SC			9-7/8" (25 cm)	15-1/3" (39 cm)	
9019-720-17SC	SCANLAN® Scanturian® MIS Super Cut™ Scissors 15° Curved Blades		6-3/4" (17 cm)	12-1/4" (31 cm)	
9019-720-25SC			9-7/8" (25 cm)	15-1/3" (39 cm)	
9019-722-17SC	SCANLAN® Scanturian® MIS Super Cut™ Scissors 30° Curved Blades		6-3/4" (17 cm)	12-1/4" (31 cm)	
9019-722-25SC			9-7/8" (25 cm)	15-1/3" (39 cm)	
9019-724-17SC	SCANLAN® Scanturian® MIS Super Cut™ Scissors 70° Curved Blades		6-3/4" (17 cm)	12-1/4" (31 cm)	
9019-724-25SC			9-7/8" (25 cm)	15-1/3" (39 cm)	






# SCANLAN® Scanturian® MIS Scissors



9019-740-17SC  
Actual Size



9019-740-25SC  
Actual Size



## Scanturian® MIS Scissors

- Well-balanced instrument
- High-quality and durable
- Super Cut™ Scissors cut cleanly and without slipping while the handle provides tactile feedback



## SCANLAN® Scanturian® MIS Scissors

Cat. No.	Description	Distal Tip	Working Length	Total Length	Handle Color
9019-740-17SC	SCANLAN® Scanturian® MIS Super Cut™ Potts Scissors 45° Angled Fine Blades		6-3/4" (17 cm)	12-1/4" (31 cm)	
9019-740-25SC			9-7/8" (25 cm)	15-1/3" (39 cm)	
9019-748-17SC	SCANLAN® Scanturian® MIS Super Cut™ Potts Scissors 90° Angled Fine Blades		6-3/4" (17 cm)	12-1/4" (31 cm)	
9019-748-25SC			9-7/8" (25 cm)	15-1/3" (39 cm)	
9019-750-17SC	SCANLAN® Scanturian® MIS Super Cut™ Potts Scissors 125° Angled Fine Blades		6-3/4" (17 cm)	12-1/4" (31 cm)	
9019-750-25SC			9-7/8" (25 cm)	15-1/3" (39 cm)	

Please Note...  
Each Scanlan instrument is individually handcrafted to provide the surgeon with the highest quality, longest-lasting, precision instrument possible.  
Reflecting this unique process, all stated dimensions are approximate.

Please check with your Scanlan Representative for Firm Sale status. Custom modifications available upon request.  
For a complete listing of models visit... [www.scanlaninternational.com](http://www.scanlaninternational.com)





SCANLAN®  
Scanturian® Sterilization Tray



SCANLAN®  
V.I.P.™ Sterilization Trays


Metal Sterilization Trays, Aluminum, Scan-Guard™ Silicone Inserts, Detachable Lids

Cat. No.	Description	Dimensions	Capacity
2160-03	SCANLAN® Scanturian® Valued Instrument Protection (V.I.P.™) Sterilization Tray	21-3/4" x 10-1/4" x 2-1/2" (55.2 cm x 26 cm x 6.4 cm)	Up to 10 MIS Instruments

All trays are steam, Immediate Use Steam Sterilization (IUSS) and ETO sterilization compatible  
Firm Sale - May be available for evaluation



SCANLAN®  
Scanturian® Sterilization Tray



SCANLAN®  
V.I.P.™ Sterilization Trays

Metal Sterilization Trays, Aluminum, Scan-Guard™ Silicone Inserts, Detachable Lids

Cat. No.	Description	Dimensions	Capacity
2161-03	SCANLAN® Scanturian® Valued Instrument Protection (V.I.P.™) Sterilization Tray with additional tray insert	21-3/4" x 10-1/4" x 4" (55.2 cm x 26 cm x 10 cm)	Up to 18 MIS Instruments

All trays are steam, Immediate Use Steam Sterilization (IUSS) and ETO sterilization compatible  
Firm Sale - May be available for evaluation







# SCANLAN® Scanturian® Sterilization Tray



SCANLAN®  
**V.I.P.™ Cleaning Tray**  
Metal Cleaning Tray, Aluminum, Scan-Guard™ Silicone Inserts

Cat. No.	Description	Dimensions	Capacity
2162-03	SCANLAN® Valued Instrument Protection (V.I.P.™) Cleaning Tray	18-7/8" x 9-7/8" x 2-1/4" (48 cm x 25 cm x 5.8 cm)	Up to 10 instruments

Firm Sale - May be available for evaluation

**For more information or to place an order,**  
please contact your local Scanlan representative, or call...

U.S. & Canada: 800-328-9458 International: 651-298-0997  
FAX: 651-298-0018  
[www.scanlaninternational.com](http://www.scanlaninternational.com)



MK-B-SCANTURIAN Rev. C